



# Coagulation Analyzer

## Micro fluidic Coagulation analysis system (Electro chemical Method)

### Fast:

Result delivered in 7 minutes

### Easy Operation:

Single drop of whole blood, One-step rapid test

### Accurate Result:

Comparable diagnostic accuracy

### Multiple parameters

✓ PT      ✓ APTT

✓ FIB      ✓ TT

\*ACT Coming soon



## Coagulation dysfunction screening



- Skin Petechiae, Purpura or Ecchymosis
- Gum bleeding
- Wound bleeding
- Gastrointestinal bleeding
- Mucosal bleeding
- Hematuria and Bloody stool
- Anemia

## Pre-operative examination



- Patients with underlying Internal disease
- Pyometra
- Tumor
- Surgery
- Biopsy

## Diagnosis and Prognosis of Internal Disease



- Multiple organ dysfunction and failure
- Obstetric diseases and complications
- Leukemia, sepsis and other blood diseases
- Immune system disease

# Operating procedures

- 1 Click test
- 2 Scan QR code
- 3 Insert reagent strip
- 4 Add sample
- 5 Start detection

